

# Capp Rondo

3 D SHAKER

**Capp Rondo 3D Shaker** operates at a variable speed range from 10 to 120 RPM. The brushless DC motor ensures stable, quiet and maintenance free operation, which is perfectly suited for extended runs. The three dimensional rocking and tumbling motion provides thorough homogenous mixing.

The user can save up to 9 individual protocols or use the pulse mode with Run-Pause-Run operation at a programmable time interval ranging from 30 to 99 seconds.

**Capp Rondo 3D Shaker** can be programmed to run from 1 up to 999 minutes, or continuously in the infinite mode.

To assure total stability on the bench, **Capp Rondo 3D Shaker** is equipped with non-slip rubber feet and comes with a spiked rubber mat. Optionally, a sticky pad mat can be ordered.

**Capp Rondo 3D Shaker** is easy to clean and complies with Protection Class IP 21.



## Specifications:

|                   |   |
|-------------------|---|
| Weight (Capacity) | 5kg   |
| Speed Range       | 10-120 RPM  |
| Motor Type        | BLDC (Brushless DC motor)                               |
| Movement          | Tumbling 7°   |
| Run Time          | 1-999 minutes or infinite                               |
| Programmable Mode | Up to 9 user defined protocols                          |
| Pulse Mode        | Run-Pause-Run programmable from 30 to 99 secs.          |
| Display           | 7 Segment LED   |
| Accessories       | Spiked rubber mat (included), Sticky pad mat (optional) |
| Plate size        | 350x350mm   |
| Protection Class  | IP 21   |

## Ordering Information

| Cat No. | Description                      |
|---------|----------------------------------|
| CR3-08X | CappRondo 3D Shaker 5kg capacity |

[www.capp.dk](http://www.capp.dk) – +45 6613 6140 – [info@capp.dk](mailto:info@capp.dk)