

 CAPP APS  
LANDBRUGSVEJ 10  
DK-5260 ODENSE  
DENMARK

TEL. (+45) 6613 6140  
FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK  
WEB: WWW.CAPP.DK

CAPP  
D E N M A R K

## Capp *Rondo* MICROCENTRIFUGE

The CAPP Rondo microcentrifuges have a lot more to offer compared to the basic microcentrifuge available in the market. It is available in two options, where the #CR-68 is the basic version with fixed speed of 6000RPM/2000G, where #CR-68X is the advanced version with programmable speed of up to 6000RPM/2000G and a timer function, all operated by two buttons and a digital display. Additional features and benefits include:

- Additional PCR strip rotor (2X8X0,2ML) and reduction adaptors for 0.2/0.4ML microtubes supplied along with product.
- Supplied with 8 slot closed rotor for quiet & powerful operation.
- Electronic safety brake provides an immediate stop on lid opening.
- Quick operation by simply closing and opening the lid.
  - Imbalance cut-off safety feature (CR-68X).
  - Digital calibration function (CR-68X).
- Digital display of speed and timer function (CR-68X).



## Ordering **Info** Capp *Rondo*

Cat No.	Description
CR-68	Capp Rondo Microcentrifuge 6000rpm/2000g
CR-68X	Capp Rondo Microcentrifuge w/ adjustable speed and timer function, max. 6000rpm/2000g